


RONGEURS

GEN2 MEDICAL Rongeurs are manufactured in Germany by instrument makers with years of experience, utilizing the highest quality stainless steel and processes. GEN2 offers one of the most comprehensive lines of rongeurs anywhere, all backed by our warranty and our commitment to quality.

We offer many different tip angles, cup sizes and lengths for your individual preference.



Rongeurs

Cat#	Description
G2N2701	ADSON RONG STR 8mm 8"
G2N2704	BACON RONG ANG 5mm 8"
G2N2705	BEYER RONG CVD 3mm 7"
G2N2707	LEMPERT RONG 2mm STR 6 1/4"
G2N2708	LEMPERT RONG 2mm CVD 6 1/4"
G2N2709	LEMPERT RONG 2mm STR 7 1/2"
G2N2710	LEMPERT RONG 2mm CVD 7 1/2"
G2N2711	RUSKIN RONG 7 1/4" STR 5mm
G2N2712	RUSKIN RONG 7 1/4" CVD 5mm
G2N2713	RUSKIN RONG 7 1/4" STR 6mm
G2N2714	RUSKIN RONG 7 1/4" CVD 6mm
G2N2715	RUSKIN RONG 6" STR 2mm
G2N2716	RUSKIN RONG 6" CVD 2mm
G2N2717	RUSKIN RONG 6" STRONG CVD 3mm
G2N2718	RUSKIN RONG 6" CVD 4mm
G2N2719	STOOKEY RONG ANG 3mm
G2N2720	STILLE-HORSLEY BONE CUTTING FCP
G2N2721	STILLE-LISTON BONE CUTTING FCP
G2N2722	STILLE-LUER RONG STR 8 1/2"
G2N2723	STILLE-LUER RONG CVD 8 1/2"
G2N2724	STILLE-LUER RONG STR 10 1/2"
G2N2725	STILLE-LUER RONG CVD 10 1/2"

Cat#	Description
G2N2726	STILLE-LUER RONG DBILL ANG 5mm 9"
G2N2727	LEKSELL RONG REG CVD 9" 8mm
G2N2728	LEKSELL RONG REG CVD 9" 5mm
G2N2729	LEKSELL RONG REG CVD 9" 3mm
G2N2730	LEKSELL RONG FULL CVD 9" 8mm
G2N2731	LEKSELL RONG FULL CVD 9" 5mm
G2N2732	LEKSELL RONG FULL CVD 9" 3mm
G2N2733	CUSHING RONG 6mm STR
G2N2734	CUSHING RONG 6mm CVD
G2N2735	LUER RONG 6mm STR 7"
G2N2736	LUER RONG 6mm CVD 7"
G2N2737	BANE RONG 4mm 7"
G2N2738	BANE RONG 5mm 7"
G2N2739	ZAUFEL-JANSEN RONGEUR 4mm 7"
G2N2740	RUSKIN RONG DBILL 5mm 7 1/4"
G2N2741	SMITH-P RONG 3mm 9 1/2" STR
G2N2742	SMITH-P RONG 3mm 9 1/2" SL-CVD
G2N2743	SMITH-P RONG 3mm 9 1/2" FL-CVD
G2N4744	ECHLIN RONG ANG 8 3/4" 2mm
G2N2745	ECHLIN RONG ANG 8 3/4" 3mm
G2N2746	ECHLIN RONG ANG 8 3/4" 4mm
G2N2747	STILLE-LUER-ECHLIN RONG ANG 9" 5mm

**MANY OTHER STYLES AND SIZES AVAILABLE

Our Rongeurs should provide trouble free use when cared for properly. All rongeurs should be inspected prior to every use, kept sharp and lubricated regularly. Special care should be taken to inspect for dullness and for any loose screws that could back out during surgery. All damaged instruments should be removed from service. All GEN2 products are guaranteed to be free of defects in workmanship and materials when used properly for their intended surgical purpose. Any product proving to be defective will be replaced or repaired, at our discretion, at no charge to the customer.